

## Laser treatment at Hands-On-Care

To facilitate healing of shoulder, elbow, wrist and hand, we at Hands-On-Care now offer Laser Treatments.

The **904 Laser** sends light (photons) to the injured tissues and stimulates healing and reduces pain instantly. No heat is generated thus:

<ul> <li>It is safe for use anywhere on the body except eyes.</li> <li>It repairs and strengthens at a remarkably fast rate.</li> <li>The treatment does NOT wear off like most medicines.</li> </ul>	Treatments last for 5-10 minutes / extremity based on the severity of the condition/ symptoms. It may take up to 6- 9 treatment sessions.
<ul> <li>It increases your body's own healing power and enhances that process.</li> <li>If you don't reinjure the area, the pain should not return.</li> <li>The LaserStim emitter has an interactive neurostimulation device that helps identify probable areas of treatments via TARGET Technology</li> </ul>	The 904 super pulsed class IV Laser treatment not only helps injured tissues and wounds heal faster, it also makes them <i>stronger and less vulnerable</i> to re-injury.
combine the entire ra	he Multi Radiance Technology <sup>™</sup> which is the first of its kind to nge of the therapeutic light spectrum and deliver therapeutic c manner. The syncing of the therapeutic wavelengths creates



Hands-On Care uses the Multi Radiance Technology<sup>m</sup> which is the first of its kind to combine the entire range of the therapeutic light spectrum and deliver therapeutic energy in a synergistic manner. The syncing of the therapeutic wavelengths creates an optimal environment that reinforces the effects of the individual wavelengths to create the world's most advanced pain relieving modality.

This **FDA** cleared the device to provide an effective non-surgical, drug-free treatment option. It features a unique combination of radiances that work synergistically for optimal pain relief and healing. The MR4 includes the patented **LaserSweep** technology which creates a constantly changing frequency to scan the selected target tissue from superficial to deep. This prevents biological adaptation to treatment for maximum impact. The MR4 excels where most laser devices fall

short. It delivers higher peak power yet is safe without the risk of thermal damage found with class IV lasers.

Synergistic Effect of MRM therapy provides benefits of light **therapy and neurostimulation**, **identifying optimal treatment sites**, and provides **highly targeted treatment** where photons need to be delivered and electro-stimulation needs to be applied.

So for faster wound healing, reducing pain and inflammation from tendonitis, or other injuries, try laser for your patients for a faster recovery.

Requires approval of workers comp. Private and Medicare patients can buy the package.

499 BLOSSOM HILL RD, SAN JOSE, CA 95120

536 408-268-8727 www

27 www.HandsOnCareTherapy.com





Results